

# Technical Data


## CIELO™ 2 BTE

### Premium Features

- Programmable BTE instruments with optional ePocket™
- For mild to moderate hearing losses
- Data logging
- High performance, automatic and adaptive directional microphone
- Automatic and adaptive feedback cancellation
- Automatic situation detection including music detection
- Adaptive noise reduction and speech enhancement in 6 channels
- Wind noise reduction
- 4 individual hearing programs (only with ePocket™) for microphone, audio shoe and/or telecoil

### Additional Features

- Alerting tones for low battery
- AutoPhone™ switchless telecoil
- Audio Input
- Volume control
- Battery compartment with lock and On/Off switch
- Battery type 13
- Nanocoated housing



### Options & Accessories

- Available in the following colors: beige, brown, grey, granite, silver, black, pearl white, light pink, light blue, transparent and translucent fun colors: purple, green, blue, orange, and pink
- Small earhook
- Audio shoe
- Materials and accessories to support pediatric fittings

## Technical Data

[www.usa.siemens.com/hearing](http://www.usa.siemens.com/hearing)

**SIEMENS**

# CIELO 2 BTE Technical Data

		2cc coupler Standard ANSI S3.22-2003	
		Standard Earhook	LifeTube*
<b>Output Sound Pressure Level</b>	Peak	123 dB	123 dB
	HF-Average OSPL 90	122 dB	114 dB
<b>Gain (Input 50 dB)</b>	Peak	55 dB	55 dB
	HF-Average	52 dB	43 dB
	Reference test gain	45 dB	37 dB
<b>Frequency Range</b>	Low frequency limit	100 Hz	100 Hz
	High frequency limit	6100 Hz	6700 Hz
<b>Total Harmonic Distortion</b>	500 Hz	4%	4%
	800 Hz	3%	3%
	1600 Hz	1%	1%
<b>Equivalent Input Noise</b>		18 dB	18 dB
<b>Inductive Coil Sensitivity</b>	HFA SPLITS** (left/right)	108/105 db	99/96 dB
	STS*** (left/right)	3/0 dB	2/-1 dB
<b>AGC-O (CK = -21 dB)</b>	Attack time	4 ms	4 ms
	Release time	100 ms	100 ms
<b>Battery</b>	Type	13 cell zinc-air	13 cell zinc-air
	Current Drain	1.2 mA	1.2 mA
	Typical Life	~190 h	~190 h

Technical information for wireless function:  
 Operating frequencies:  $f_{\text{low}} = 115 \text{ kHz}$ ,  $f_{\text{high}} = 120 \text{ kHz}$ ;  
 Rated H-field strength (maximum):  $-11,5 \mu\text{A/m}$  at 3 meters.

\* Measured with S-LifeTube, Size 2

\*\*SPLITS = Coupler SPL for an Inductive Telephone Simulator


\*\*\*STS = Simulated Telephone Sensitivity

Measure instructions: Instrument in linear setting.


# CIELO 2 BTE Basic Data

## Basic Data w/Standard Earhook


Output Sound Pressure Level ( $L_i = 90$  dB)  
 Maximum Gain ( $L_i = 50$  dB)




Frequency Response ( $L_i = 60$  dB)  
 Basic Frequency Response ( $L_i = 60$  dB)



Effect of MPO (FOG,  $f = 2$  kHz)




Inductive Response Left ( $H = 31.6$  mA/m)  
 Inductive Response Right ( $H = 31.6$  mA/m)




## Basic Data with LifeTube™


Output Sound Pressure Level ( $L_i = 90$  dB)  
 Maximum Gain ( $L_i = 50$  dB)




Frequency Response ( $L_i = 60$  dB)  
 Basic Frequency Response ( $L_i = 60$  dB)










Effect of MPO (FOG,  $f = 2$  kHz)



Inductive Response Left ( $H = 31.6$  mA/m)  
 Inductive Response Right ( $H = 31.6$  mA/m)



# CIELO 2 BTE Fitting Information

	<h3><u>Programming Socket</u></h3> <p>The programming socket lies under a flap below the volume control. Use a suitable tool to open the flap. After the programming procedure is complete, close the flap using your fingernail.</p>
	<h3><u>Changing the Typeplates</u></h3> <p>To mark the left and right side of the CIELO 2 instruments when fitting binaurally, exchange the housing colored typeplates on the inner curve of the instruments for blue (left) and red (right). Use a suitable tool to lift and remove the typeplate. Lock the two pins on the new typeplate into the openings and press gently into position with your finger.</p>
	<h3><u>ePocket™</u></h3> <p>The CIELO 2 hearing system supports the use of ePocket™, a bi-directional remote control with read out function. ePocket can change the hearing programs and the volume of CIELO 2. The ePocket read out function will display the current program, volume level and battery status of the instrument(s). ePocket includes a cover and clip.</p>
	<h3><u>Audio Shoe</u></h3> <p>To attach an audio shoe to CIELO 2, open the battery compartment to the first stop. Open the flap under the programming socket and hook the audio shoe in from the front so the gold contacts on the audio shoe are aligned with the gold contacts under the instrument flap. Press the audio shoe into place so that it attaches to the back and the curve of the audio shoe presses against the bottom of the hearing instrument. To remove the audio shoe, pull it toward the front. Then unhook the audio shoe and close the flap.</p>
	<h3><u>Earhook, small</u></h3> <p>For optimal fitting to smaller ears, a small earhook is an available option.</p>
	<h3><u>Open Fitting Option</u></h3> <p>Slim S-LifeTubes and soft, flexible LifeTips are available in various sizes to conduct an open fitting for select, Siemens CIELO 2 model BTEs.</p>
	<h3><u>Pediatric Accessories</u></h3> <p>A Use and Care kit containing items to care for hearing instruments (blower, stethoset, drying kit, and more) and a teddy bear wearing hearing instruments packed into a colorful lunchbox; storybooks for children about hearing loss, and other materials are available to support pediatric fittings. For more information or to order pediatric materials, contact your Siemens Sales Representative.</p>

## Siemens Hearing Instruments, Inc. locations

### United States Headquarters/Northeast Manufacturing Facility:

10 Constitution Avenue, P.O. Box 1397, Piscataway, NJ 08855-1397 • (732) 562-6600 or (800) 766-4500

**Midwest Sales Facility:** (847) 808-1200 or (800) 333-9083

**West Manufacturing Facility:** (562) 404-4531 or (800) 998-9787

**Technical Support for Software and Systems:** (888) 231-1333

[www.usa.siemens.com/hearing](http://www.usa.siemens.com/hearing)

### Siemens Hearing Instruments A Division of Siemens Canada Limited

320 Pinebush Road, Cambridge, Ontario, Canada N1T 1Z6 • (519) 622-5200 or (800) 663-0620

[www.siemens.ca/hearing](http://www.siemens.ca/hearing)

Information in this brochure subject to change without notice. © 2007 Siemens Hearing Instruments, Inc. All rights reserved