

SIEMENS Acuris

Acuris – the world's most advanced hearing aid system



Now, for the first time, wearing two hearing instruments offers a new level of sophistication not possible before. Because of E2E (ear-to-ear) wireless technology, *ACURIS* hearing instruments actually communicate with each other. They sense, analyze, interpret, share and adjust continuously and automatically.



Siemens' goal in developing *ACURIS* was to ensure the end result produced incredible performance with a new level of ease, simplicity and comfort. For most people that means maximizing the real power and benefit of wearing two instruments – performing as one.



- Press the listening mode button or turn the volume control on one hearing aid and both sides will change simultaneously.
- No more adjusting or balancing one hearing instrument, then the other, to keep both in sync.
- The special *ACURIS* ePocket remote control is an option that both adjusts and displays the hearing instrument's settings in your hand.



Product Highlights:

- E2E wireless functions – the instruments work in total harmony together
- Extremely advanced Digital Signal Processing with 16 bands/8 channels
- High-speed Feedback Cancellation...very effective
- Remote control option for maximum control of even the tiniest models



Custom in-the-ear models

Full-shell	\$2,000
Half-shell	\$2,000
Canal/Mini Canal	\$2,000
C.I.C.	\$2,100

BTE models

Acuris S (Std.)	\$2,000
Acuris P (Hi-Power)	\$2,100
Optional ePocket Remote....	\$150

Additional options

(for custom I-T-E models)

*Telephone coil.....	\$50
Switchless T-coil.....	\$90
*Directional mics	\$200

*Standard and no charge on BTE models